

Initiating Search

November 10, 2025, 3:54 PM

• Search:

CAS Lexicon Search:

N-Chlorosuccinimide - Preferred Concept

Search Tasks

Task	Result Type	View
Returned Reference Results from CAS Lexicon (4,885)	 References	View Results
Exported: Retrieved Related Reaction Results + Filters (2,534)	 Reactions	View Results
Filtered By:		
Yield:	90-100%, 80-89%, 70-79%	
Document Type:	Journal	

Copyright © 2025 American Chemical Society (ACS). All Rights Reserved.


Internal use only. Redistribution is subject to the terms of your CAS SciFinder License Agreement and CAS information Use Policies.

Reactions (50)

[View in CAS SciFinder](#)

Scheme 1 (1 Reaction)

Steps: 1 Yield: 100%


 Suppliers (99)

 Suppliers (75)

 Suppliers (5)

31-084-CAS-8970901

Steps: 1 Yield: 100%

Bromination and chlorination of pyrrole and some reactive 1-substituted pyrroles


By: Gilow, H. M.; et al

Journal of Organic Chemistry (1981), 46(11), 2221-5.

No Data Available

Scheme 2 (1 Reaction)

Steps: 1 Yield: 100%


 Suppliers (59)

 Suppliers (81)

31-084-CAS-3949451

Steps: 1 Yield: 100%

Effect of N-halosuccinimides on 5-, 7- and 5,7-8-quinolinol sulfonic acids


By: Gershon, Herman; et al

Trends in Heterocyclic Chemistry (2007), 12, 67-73.

1.1 **Reagents:** N-Bromosuccinimide
Solvents: Acetic acid, Acetonitrile; 10 min, 60 °C

Scheme 3 (1 Reaction)

Steps: 1 Yield: 100%


 Suppliers (4)

 Supplier (1)

31-084-CAS-5539733

Steps: 1 Yield: 100%

Reactions of N-aminopyrazoles with halogenating reagents and synthesis of 1,2,3-triazines

By: Osawa, Akio; et al

Chemical & Pharmaceutical Bulletin (1988), 36(10), 3838-48.

1.1 **Reagents:** Bromine
Solvents: Dichloromethane

Scheme 4 (1 Reaction)

Steps: 1 Yield: 100%



31-084-CAS-13866343

Steps: 1 Yield: 100%


Regioselective approach to multi-substituted benzenes

1.1 Reagents: *N*-Bromosuccinimide
 Solvents: Dichloromethane; 2.5 h, 0 °C

By: Seo, Hana; et al
 Chemistry Letters (2011), 40(7), 744-746.

Scheme 5 (1 Reaction)

Steps: 1 Yield: 100%



Supplier (1)

31-084-CAS-2821953

Steps: 1 Yield: 100%


Design and synthesis of new tricyclic indoles as potent modulators of the S1P₁ receptor

1.1 Reagents: Chlorosuccinimide
 Solvents: Dichloromethane; 0 °C → rt

By: Buzard, Daniel J.; et al
 Bioorganic & Medicinal Chemistry Letters (2015), 25(3), 659-663.

Scheme 6 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (86)

Suppliers (100)

31-084-CAS-23031761

Steps: 1 Yield: 100%

Activator free, expeditious and eco-friendly chlorination of activated arenes by *N*-chloro-*N*-(phenylsulfonyl)benzene sulfonamide (NCBSI)

1.1 Reagents: *N*-Chloro-*N*-(phenylsulfonyl)benzenesulfonamide
 Solvents: Acetonitrile; 5 min, 20 - 25 °C


1.2 Reagents: Water
 Solvents: Dichloromethane; 10 - 15 min, 20 - 25 °C

By: Misal, Balu; et al
 Tetrahedron Letters (2021), 63, 152689.


Experimental Protocols

Scheme 7 (1 Reaction)


Steps: 1 Yield: 100%




31-084-CAS-15121543 Steps: 1 Yield: 100%	Consecutive Thiophene-Annulation Approach to π-Extended Thienoacene-Based Organic Semiconductors with [1]Benzot hieno[3,2-b][1]benzothiophene (BTBT) Substructure
1.1 Reagents: Butyllithium Solvents: Tetrahydrofuran, Hexane; 0 °C; 0 °C → rt; 1 h, rt 1.2 Reagents: 1,2-Dibromo-1,1,2,2-tetrachloroethane; 0 °C; 6 h, rt	By: Mori, Takamichi; et al
Experimental Protocols	Journal of the American Chemical Society (2013), 135(37), 13900-13913.

Scheme 8 (1 Reaction)	Steps: 1 Yield: 100%
 Suppliers (99) Suppliers (87)	

31-084-CAS-23033350 Steps: 1 Yield: 100%	Activator free, expeditious and eco-friendly chlorination of activated arenes by N-chloro-N-(phenylsulfonyl)benzene sulfonamide (NCBSI)
1.1 Reagents: <i>N</i> -Chloro- <i>N</i> -(phenylsulfonyl)benzenesulfonamide Solvents: Acetonitrile; 10 min, 20 - 25 °C 1.2 Reagents: Water Solvents: Dichloromethane; 10 - 15 min, 20 - 25 °C	By: Misal, Balu; et al Tetrahedron Letters (2021), 63, 152689.
Experimental Protocols	

Scheme 9 (1 Reaction)	Steps: 1 Yield: 100%
 Suppliers (44) Suppliers (2)	


31-084-CAS-9837269 Steps: 1 Yield: 100%	Development of serine protease inhibitors displaying a multicentered short (<2.3 Å) hydrogen bond binding mode: Inhibitors of urokinase-type plasminogen activator and factor Xa
1.1 Reagents: <i>N</i> -Bromosuccinimide Solvents: Acetic acid; 2 h, 85 °C	By: Verner, Erik; et al
Experimental Protocols	Journal of Medicinal Chemistry (2001), 44(17), 2753-2771.

Scheme 10 (1 Reaction)	Steps: 1 Yield: 100%
 Supplier (1)	

31-084-CAS-17240926 Steps: 1 Yield: 100% 1.1 Reagents: <i>N</i>-Bromosuccinimide Catalysts: Palladium diacetate Solvents: Chlorobenzene; 24 h, 110 °C Experimental Protocols	Palladium-Catalyzed Directed Halogenation of Bipyridine N-Oxides By: Zucker, Sina P.; et al Journal of Organic Chemistry (2017), 82(11), 5616-5635.
--	--

Scheme 11 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (69)


Suppliers (101)

Supplier (1)

31-084-CAS-10449571 No Data Available	Steps: 1 Yield: 100% Bromination and chlorination of pyrrole and some reactive 1-substituted pyrroles By: Gilow, H. M.; et al Journal of Organic Chemistry (1981), 46(11), 2221-5.
---	--

Scheme 12 (1 Reaction)

Steps: 1 Yield: 100%




Supplier (1)

31-084-CAS-8935179 Steps: 1 Yield: 100% 1.1 Reagents: Chlorosuccinimide Solvents: Dichloromethane; 0 °C → rt	Design and synthesis of new tricyclic indoles as potent modulators of the S1P₁ receptor By: Buzard, Daniel J.; et al Bioorganic & Medicinal Chemistry Letters (2015), 25(3), 659-663.
---	---

Scheme 13 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (99)

31-084-CAS-18183871 Steps: 1 Yield: 100% 1.1 Solvents: Dichloromethane; 0 °C → rt	Design and synthesis of new tricyclic indoles as potent modulators of the 51P₁ receptor By: Buzard, Daniel J.; et al Bioorganic & Medicinal Chemistry Letters (2015), 25(3), 659-663.
---	--

Scheme 14 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (77)

31-084-CAS-19489252 Steps: 1 Yield: 100% 1.1 Reagents: Potassium bromide, Vanadate (VO₄³⁻), sodium (1:3), (7-4)- Catalysts: Flavin mononucleotide Solvents: Acetonitrile, Water; 16 h, pH 6, 25 °C Experimental Protocols	Atom-Economic Electron Donors for Photobiocatalytic Halogenations By: Seel, Catharina Julia; et al ChemCatChem (2018), 10(18), 3960-3963.
--	---

Scheme 15 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (29)

31-084-CAS-22871968 Steps: 1 Yield: 100% 1.1 Reagents: Acetic acid, <i>N</i>-Iodosuccinimide Catalysts: Silver triflate, Iridium, di-μ-chlorodichlorobis[(1,2,3,4,5-η)-1,2,3,4,5-pentamethyl-2,4-cyclopentadien-1-yl]di- Solvents: 1,2-Dichloroethane; 1 h, 60 °C Experimental Protocols	Identification of key functionalization species in the Cp*Ir(III)-catalyzed-ortho halogenation of benzamides By: Guzman Santiago, Alexis J.; et al Dalton Transactions (2020), 49(45), 16166-16174.
---	---

Scheme 16 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (99)

31-084-CAS-18177038	Steps: 1 Yield: 100%	Design and synthesis of new tricyclic indoles as potent modulators of the 51P₁ receptor By: Buzard, Daniel J.; et al Bioorganic & Medicinal Chemistry Letters (2015), 25(3), 659-663.
1.1 Solvents: Dichloromethane; 0 °C → rt		

Scheme 17 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (88)

Suppliers (60)

31-084-CAS-8621113	Steps: 1 Yield: 100%	Efficient Oxidative Chlorination of Aromatics on Saturated Sodium Chloride Solution By: Gu, Liuqun; et al Advanced Synthesis & Catalysis (2013), 355(6), 1077-1082.
1.1 Reagents: Sodium chloride, Potassium persulfate Solvents: Acetonitrile, Water; 1.5 h, 100 °C		
1.2 Reagents: Sodium thiosulfate Solvents: Water		
Experimental Protocols		

Scheme 18 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (71)


Suppliers (101)

Suppliers (2)

31-084-CAS-7983078	Steps: 1 Yield: 100%	Bromination and chlorination of pyrrole and some reactive 1-substituted pyrroles By: Gilow, H. M.; et al Journal of Organic Chemistry (1981), 46(11), 2221-5.
No Data Available		

Scheme 19 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (62)

Suppliers (41)

31-614-CAS-39834437 Steps: 1 Yield: 100%	Exploring boron applications in modern agriculture: A structure-activity relationship study of a novel series of multi-substitution benzoxaboroles for identification of potential fungicides By: Liu, Chunliang; et al Bioorganic & Medicinal Chemistry Letters (2021), 43, 128089.
1.1 Reagents: Chlorosuccinimide Solvents: Dimethylformamide; 3 h, rt 1.2 Reagents: Water; 2 min, 0 °C	
Experimental Protocols	

Scheme 20 (1 Reaction)


Steps: 1 Yield: 100%



31-084-CAS-20103702 Steps: 1 Yield: 100%	Aromatic Halogenation Using N-Halosuccinimide and PhSSi Me₃ or PhSSPh By: Hirose, Yuuka; et al Journal of Organic Chemistry (2019), 84(11), 7405-7410.
1.1 Reagents: <i>N</i> -Bromosuccinimide Catalysts: Diphenyl disulfide Solvents: Acetonitrile; 12 h, rt 1.2 Reagents: Sodium bicarbonate, Dipotassium phosphate Solvents: Water; rt	
Experimental Protocols	

Scheme 21 (1 Reaction)


Steps: 1 Yield: 100%



31-084-CAS-19489247 Steps: 1 Yield: 100%	Atom-Economic Electron Donors for Photobiocatalytic Halogenations By: Seel, Catharina Julia; et al ChemCatChem (2018), 10(18), 3960-3963.
1.1 Reagents: Potassium bromide, Vanadate (VO ₄ ³⁻), sodium (1:3), (7-4)- Catalysts: Flavin mononucleotide, Haloperoxidase Solvents: Acetonitrile, Water; 16 h, pH 6, 25 °C	
Experimental Protocols	

Scheme 22 (1 Reaction)


Steps: 1 Yield: 100%



31-084-CAS-23032184 Steps: 1 Yield: 100%	Activator free, expeditious and eco-friendly chlorination of activated arenes by N-chloro-N-(phenylsulfonyl)benzene sulfonamide (NCBSI) By: Misal, Balu; et al Tetrahedron Letters (2021), 63, 152689.
1.1 Reagents: <i>N</i> -Chloro- <i>N</i> -(phenylsulfonyl)benzenesulfonamide Solvents: Acetonitrile; 20 min, 20 - 25 °C	
1.2 Reagents: Water Solvents: Dichloromethane; 10 - 15 min, 20 - 25 °C	
Experimental Protocols	

Scheme 23 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (133)

Suppliers (86)

31-084-CAS-23029533 Steps: 1 Yield: 100%	Activator free, expeditious and eco-friendly chlorination of activated arenes by N-chloro-N-(phenylsulfonyl)benzene sulfonamide (NCBSI) By: Misal, Balu; et al Tetrahedron Letters (2021), 63, 152689.
1.1 Reagents: <i>N</i> -Chloro- <i>N</i> -(phenylsulfonyl)benzenesulfonamide Solvents: Acetonitrile; 10 min, 20 - 25 °C	
1.2 Reagents: Water Solvents: Dichloromethane; 10 - 15 min, 20 - 25 °C	
Experimental Protocols	

Scheme 24 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (61)

Suppliers (76)

31-084-CAS-6497966 Steps: 1 Yield: 100%	Efficient Oxidative Chlorination of Aromatics on Saturated Sodium Chloride Solution By: Gu, Liuqun; et al Advanced Synthesis & Catalysis (2013), 355(6), 1077-1082.
1.1 Reagents: Sodium chloride, Potassium persulfate Solvents: Acetonitrile, Water; 1.5 h, 100 °C	
1.2 Reagents: Sodium thiosulfate Solvents: Water	
Experimental Protocols	

Scheme 25 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (113)

Suppliers (48)

31-084-CAS-6770321	Steps: 1 Yield: 100%	Step-economic synthesis of (±)-debromoflustramine A using indole C3 activation strategy
1.1 Reagents: 1,4-Dimethylpiperazine, Chlorosuccinimide Solvents: Dichloromethane; 2 h, 0 °C		By: Ignatenko, Vasily A.; et al
1.2 Reagents: Trichloroacetic acid, 3-Methyl-2-buten-1-ol Solvents: Dichloromethane; 1 h, 0 °C; overnight, rt		Tetrahedron Letters (2011), 52(12), 1269-1272.
Experimental Protocols		

Scheme 26 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (85)

Suppliers (54)

31-084-CAS-19489243	Steps: 1 Yield: 100%	Atom-Economic Electron Donors for Photobiocatalytic Halogenations
1.1 Reagents: Potassium bromide, Vanadate (VO ₄ ³⁻), sodium (1:3), (7-4)- Catalysts: Flavin mononucleotide, Haloperoxidase Solvents: Acetonitrile, Water; 16 h, pH 6, 25 °C		By: Seel, Catharina Julia; et al
Experimental Protocols		ChemCatChem (2018), 10(18), 3960-3963.

Scheme 27 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (2)

31-614-CAS-24202727	Steps: 1 Yield: 100%	Design of novel, water soluble and highly luminescent europium labels with potential to enhance immunoassay sensitivities
1.1 Reagents: Chlorosuccinimide, Sodium iodide Solvents: Acetic acid; 1.5 h, rt		By: Sund, Henri; et al
1.2 Reagents: Sodium bicarbonate Solvents: Water; neutralized, rt		Molecules (2017), 22(10), 1807/1-1807/23.
Experimental Protocols		

Scheme 28 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (88)

Suppliers (58)

31-084-CAS-20103708	Steps: 1 Yield: 100%	Aromatic Halogenation Using N-Halosuccinimide and PhSSi Me₃ or PhSSPh By: Hirose, Yuuka; et al Journal of Organic Chemistry (2019), 84(11), 7405-7410.
1.1 Reagents: N-Bromosuccinimide Catalysts: Diphenyl disulfide Solvents: Acetonitrile; 1 h, rt		
1.2 Reagents: Sodium bicarbonate, Dipotassium phosphate Solvents: Water; rt		

Scheme 29 (1 Reaction)


Steps: 1 Yield: 100%



31-084-CAS-19489250	Steps: 1 Yield: 100%	Atom-Economic Electron Donors for Photobiocatalytic Halogenations By: Seel, Catharina Julia; et al ChemCatChem (2018), 10(18), 3960-3963.
1.1 Reagents: Potassium bromide, Vanadate (VO ₄ ³⁻), sodium (1:3), (7-4)- Catalysts: Flavin mononucleotide Solvents: Acetonitrile, Water; 16 h, pH 6, 25 °C		
Experimental Protocols		

Scheme 30 (1 Reaction)


Steps: 1 Yield: 100%



31-084-CAS-15017849	Steps: 1 Yield: 100%	Efficient Oxidative Chlorination of Aromatics on Saturated Sodium Chloride Solution By: Gu, Liuqun; et al Advanced Synthesis & Catalysis (2013), 355(6), 1077-1082.
1.1 Reagents: Sodium chloride, Potassium persulfate Solvents: Acetonitrile, Water; 1 h, 100 °C		
1.2 Reagents: Sodium thiosulfate Solvents: Water		
Experimental Protocols		

Scheme 31 (1 Reaction)


Steps: 1 Yield: 100%



31-084-CAS-20103710 Steps: 1 Yield: 100% 1.1 Reagents: <i>N</i>-Bromosuccinimide Catalysts: Diphenyl disulfide Solvents: Acetonitrile; 1 h, rt 1.2 Reagents: Sodium bicarbonate, Dipotassium phosphate Solvents: Water; rt	Aromatic Halogenation Using <i>N</i>-Halosuccinimide and PhSSiMe₃ or PhSSPh By: Hirose, Yuuka; et al Journal of Organic Chemistry (2019), 84(11), 7405-7410.
--	---

Scheme 32 (1 Reaction)


Steps: 1 Yield: 100%



31-084-CAS-19489241 Steps: 1 Yield: 100% 1.1 Reagents: Potassium bromide, Vanadate (VO₄³⁻), sodium (1:3), (7-4)- Catalysts: Flavin mononucleotide, Haloperoxidase Solvents: Acetonitrile, Water; 6 h, pH 6, 25 °C Experimental Protocols	Atom-Economic Electron Donors for Photobiocatalytic Halogenations By: Seel, Catharina Julia; et al ChemCatChem (2018), 10(18), 3960-3963.
--	---

Scheme 33 (1 Reaction)


Steps: 1 Yield: 99%



31-084-CAS-20456593 Steps: 1 Yield: 99% 1.1 Reagents: Acetic acid, <i>N</i>-Iodosuccinimide Catalysts: Rhodium(2+), tris(acetonitrile)[(1,2,3,4,5-η)-1,2,3,4,5-pentamethyl-2,4-cyclopentadien-1-yl]-, (OC-6-11)-hexafluoroantimonate(1-) (1:2) Solvents: 1,2-Dichloroethane; -196 °C; -196 °C → rt; 18 h, 60 °C	Transition metal catalyzed C7 and ortho- selective halogenation of 2-arylbenzo[d]oxazoles By: Hong, Xi; et al Organic Chemistry Frontiers (2019), 6(13), 2226-2233.
--	---

Scheme 34 (1 Reaction)


Steps: 1 Yield: 99%



31-084-CAS-8441062	Steps: 1 Yield: 99%	Synthesis of (ω-haloperfluoroalkyl)imidazoles and their derivatives By: Fujii, Shozo; et al Nagoya Kogyo Gijutsu Shikensho Hokoku (1986), 35(3), 117-29.
1.1 Reagents: Iodine		

Scheme 35 (1 Reaction)

Steps: 1 Yield: 99%




Supplier (1)

31-084-CAS-16741027	Steps: 1 Yield: 99%	Selective C(sp²)-H Halogenation of "Click" 4-Aryl-1,2,3-triazoles By: Goitia, Asier; et al Organic Letters (2017), 19(4), 962-965.
1.1 Reagents: <i>N</i> -Bromosuccinimide Solvents: 1,2-Dichloroethane; 18 h, rt		
Experimental Protocols		

Scheme 36 (1 Reaction)

Steps: 1 Yield: 99%




Suppliers (103)

Suppliers (61)

31-084-CAS-10038012	Steps: 1 Yield: 99%	Halogenation of aromatic compounds by <i>N</i>-chloro-, <i>N</i>-bromo-, and <i>N</i>-iodosuccinimide By: Tanemura, Kiyoshi; et al Chemistry Letters (2003), 32(10), 932-933.
1.1 Reagents: <i>N</i> -Iodosuccinimide Catalysts: Iron chloride (FeCl ₃) Solvents: Acetonitrile; 0.5 h, 25 °C		

Scheme 37 (1 Reaction)

Steps: 1 Yield: 99%




Suppliers (83)

Suppliers (82)

31-084-CAS-13835539	Steps: 1 Yield: 99%	Halogenation of pyrazoles using <i>N</i>-halosuccinimides in CCl₄ and in water By: Zhao, Zhi-Gang; et al Synthetic Communications (2007), 37(1), 137-147.
1.1 Reagents: <i>N</i> -Bromosuccinimide Solvents: Carbon tetrachloride; 2 h, 25 °C		
Experimental Protocols		

Scheme 38 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (75)

Suppliers (20)

31-084-CAS-23028293

Steps: 1 Yield: 99%

Activator free, expeditious and eco-friendly chlorination of activated arenes by N-chloro-N-(phenylsulfonyl)benzene sulfonamide (NCBSI)

1.1 **Reagents:** *N*-Chloro-*N*-(phenylsulfonyl)benzenesulfonamide
Solvents: Acetonitrile; 10 min, 20 - 25 °C

1.2 **Reagents:** Water
Solvents: Dichloromethane; 10 - 15 min, 20 - 25 °C


By: Misal, Balu; et al

Tetrahedron Letters (2021), 63, 152689.

Experimental Protocols

Scheme 39 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (97)

Suppliers (81)

31-084-CAS-10073755

Steps: 1 Yield: 99%

Gold-catalyzed halogenation of aromatics by N-halosuccinimides

1.1 **Reagents:** *N*-Bromosuccinimide
Catalysts: Gold trichloride
Solvents: 1,2-Dichloroethane; 5 h, 60 °C


By: Mo, Fanyang; et al

Angewandte Chemie, International Edition (2010), 49(11), 2028-2032, S2028/1-S2028/10.

Experimental Protocols

Scheme 40 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (99)

Suppliers (93)

31-084-CAS-23029114

Steps: 1 Yield: 99%

Activator free, expeditious and eco-friendly chlorination of activated arenes by N-chloro-N-(phenylsulfonyl)benzene sulfonamide (NCBSI)

1.1 **Reagents:** *N*-Chloro-*N*-(phenylsulfonyl)benzenesulfonamide
Solvents: Acetonitrile; 10 min, 20 - 25 °C

1.2 **Reagents:** Water
Solvents: Dichloromethane; 10 - 15 min, 20 - 25 °C


By: Misal, Balu; et al

Tetrahedron Letters (2021), 63, 152689.

Experimental Protocols

Scheme 41 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (90)

Suppliers (99)

31-084-CAS-23030492

Steps: 1 Yield: 99%

Activator free, expeditious and eco-friendly chlorination of activated arenes by N-chloro-N-(phenylsulfonyl)benzene sulfonamide (NCBSI)

1.1 **Reagents:** *N*-Chloro-*N*-(phenylsulfonyl)benzenesulfonamide
Solvents: Acetonitrile; 15 min, 20 - 25 °C

1.2 **Reagents:** Water
Solvents: Dichloromethane; 10 - 15 min, 20 - 25 °C


By: Misal, Balu; et al

Tetrahedron Letters (2021), 63, 152689.

Experimental Protocols

Scheme 42 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (115)

Suppliers (74)

31-084-CAS-667497

Steps: 1 Yield: 99%

Gold-catalyzed halogenation of aromatics by N-halosuccinimides

1.1 **Reagents:** *N*-Bromosuccinimide
Catalysts: Gold trichloride
Solvents: 1,2-Dichloroethane; 30 h, rt


By: Mo, Fanyang; et al

Angewandte Chemie, International Edition (2010), 49(11), 2028-2032, S2028/1-S2028/10.

Experimental Protocols

Scheme 43 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (93)

Suppliers (63)

31-084-CAS-23027431

Steps: 1 Yield: 99%

Activator free, expeditious and eco-friendly chlorination of activated arenes by N-chloro-N-(phenylsulfonyl)benzene sulfonamide (NCBSI)

1.1 **Reagents:** *N*-Chloro-*N*-(phenylsulfonyl)benzenesulfonamide
Solvents: Acetonitrile; 5 min, 20 - 25 °C

1.2 **Reagents:** Water
Solvents: Dichloromethane; 10 - 15 min, 20 - 25 °C


By: Misal, Balu; et al

Tetrahedron Letters (2021), 63, 152689.

Experimental Protocols

Scheme 44 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (93)

Suppliers (84)

31-084-CAS-23033407

Steps: 1 Yield: 99%

Activator free, expeditious and eco-friendly chlorination of activated arenes by N-chloro-N-(phenylsulfonyl)benzenesulfonamide (NCBSI)

By: Misal, Balu; et al

Tetrahedron Letters (2021), 63, 152689.


1.1 Reagents: *N*-Chloro-*N*-(phenylsulfonyl)benzenesulfonamide
Solvents: Acetonitrile; 8 min, 20 - 25 °C

1.2 Reagents: Water
Solvents: Dichloromethane; 10 - 15 min, 20 - 25 °C

Experimental Protocols

Scheme 45 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (83)

Suppliers (37)

31-084-CAS-2567980

Steps: 1 Yield: 99%

Halogenation of pyrazoles using N-halosuccinimides in CCl₄ and in water

By: Zhao, Zhi-Gang; et al


Synthetic Communications (2007), 37(1), 137-147.

1.1 Reagents: Chlorosuccinimide
Solvents: Carbon tetrachloride; 1 h, 50 °C

Experimental Protocols

Scheme 46 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (93)

Suppliers (82)

31-084-CAS-9077272

Steps: 1 Yield: 99%

Silver(I)-Catalyzed Iodination of Arenes: Tuning the Lewis Acidity of N-Iodosuccinimide Activation

By: Racys, Daugirdas T.; et al


Journal of Organic Chemistry (2016), 81(3), 772-780.

1.1 Reagents: *N*-Iodosuccinimide
Catalysts: [1,1,1-Trifluoro-*N*-[(trifluoromethyl)sulfonyl-κO]methanesulfonamidato-κO]silver
Solvents: Dichloromethane; 1.5 h, 40 °C

Experimental Protocols

Scheme 47 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (85)

Suppliers (88)

31-084-CAS-9754035

Steps: 1 Yield: 99%

Halogenation of aromatic compounds by N-chloro-, N-bromo-, and N-iodosuccinimide


- 1.1 **Reagents:** *N*-Iodosuccinimide
Catalysts: Iron chloride (FeCl₃)
Solvents: Acetonitrile; 0.5 h, 25 °C

By: Tanemura, Kiyoshi; et al

Chemistry Letters (2003), 32(10), 932-933.

Scheme 48 (3 Reactions)

Steps: 1 Yield: 99%



Suppliers (85)

Suppliers (71)

31-614-CAS-44264822

Steps: 1 Yield: 99%

Dual activation with Fe(NTf₂)₃ and I₂: Regioselective halogenation of aromatic compounds using N-halosuccinimides (NX S)

- 1.1 **Catalysts:** Iron chloride (FeCl₃), 1-Butyl-3-methylimidazolium bis(trifluoromethylsulfonyl)imide; 30 min, rt
 1.2 **Reagents:** *N*-Bromosuccinimide
Catalysts: Iodine
Solvents: Acetonitrile; 0.5 h, rt

By: Tanemura, Kiyoshi

Tetrahedron (2025), 176, 134544.

Experimental Protocols

31-614-CAS-33508000

Steps: 1 Yield: 99%

Hypervalent Chalcogenonium... π Bonding Catalysis

- 1.1 **Reagents:** *N*-Bromosuccinimide
Catalysts: 3*H*-2,1-Benzoxaselenolium, 3,3-dimethyl-1-phenyl-, tetrafluoroborate(1-) (1:1)
Solvents: Dichloromethane; 23 °C

By: Zhang, Qingyu; et al

Angewandte Chemie, International Edition (2022), 61(35), e202208009.

Experimental Protocols

31-084-CAS-3061803

Steps: 1 Yield: 99%

Gold-catalyzed halogenation of aromatics by N-halosuccinimides

- 1.1 **Reagents:** *N*-Bromosuccinimide
Catalysts: Gold trichloride
Solvents: 1,2-Dichloroethane; 1 h, rt

By: Mo, Fanyang; et al

Angewandte Chemie, International Edition (2010), 49(11), 2028-2032, S2028/1-S2028/10.

Experimental Protocols