

Initiating Search

November 4, 2025, 11:31 AM

• Search:

CAS Lexicon Search:

Stille coupling reaction - Preferred Concept

Search Tasks

Task	Result Type	View
Returned Reference Results from CAS Lexicon (4,657)	 References	View Results
Exported: Retrieved Related Reaction Results + Filters (3,197)	 Reactions	View Results
Filtered By:		
Yield:	90-100%, 80-89%	
Language:	English	

Copyright © 2025 American Chemical Society (ACS). All Rights Reserved.


Internal use only. Redistribution is subject to the terms of your CAS SciFinder License Agreement and CAS information Use Policies.

Reactions (50)

[View in CAS SciFinder](#)

Scheme 1 (1 Reaction)

Steps: 1 Yield: 100%


 Suppliers (36)

 Suppliers (25)

31-130-CAS-9716226

Steps: 1 Yield: 100%

1.1 Catalysts: Tetrakis(triphenylphosphine)palladium
Solvents: Dimethylformamide; 2 h, 110 °C


Stille Cross-Coupling for the Functionalization of the Imidazole Backbone: Revisit, Improvement, and Applications of the Method

By: Sandtorv, Alexander H.; et al

European Journal of Organic Chemistry (2015), 2015(16), 3506-3512.

Scheme 2 (1 Reaction)

Steps: 1 Yield: 100%


 Suppliers (3)

 Suppliers (53)

31-130-CAS-15803527

Steps: 1 Yield: 100%

1.1 Solvents: Toluene; 5 min, 100 °C; 20 min, 170 °C

Experimental Protocols


Utility of a heterogeneous palladium catalyst for the synthesis of a molecular semiconductor via Stille, Suzuki, and direct heteroarylation cross-coupling reactions

By: McAfee, Seth M.; et al

RSC Advances (2015), 5(33), 26097-26106.

Scheme 3 (1 Reaction)

Steps: 1 Yield: 100%



 Suppliers (2)

 Suppliers (2)

31-130-CAS-68496 Steps: 1 Yield: 100%	A C-Ring regioisomer of the marine alkaloid meridine exhibits selective in vitro cytotoxicity for solid tumours
1.1 Catalysts: Palladium, bis(acetonitrile)dichloro- Solvents: 1,4-Dioxane	By: Angel de la Fuente, J.; et al Bioorganic & Medicinal Chemistry (2001), 9(7), 1807-1814.

Scheme 4 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-9354993 Steps: 1 Yield: 100%	Synthesis of asymmetric 3,5-diaryl-4H-1,2,6-thiadiazin-4-ones via Suzuki-Miyaura and Stille coupling reactions
1.1 Catalysts: Dichlorobis(triphenylphosphine)palladium Solvents: Acetonitrile; 20 °C; 3.8 h, 82 °C	By: Ioannidou, Heraklidia A.; et al Tetrahedron (2012), 68(36), 7380-7385.

Scheme 5 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-14308511 Steps: 1 Yield: 100%	Stille Cross-Coupling for the Functionalization of the Imidazole Backbone: Revisit, Improvement, and Applications of the Method
1.1 Catalysts: Cuprous iodide, Tetrakis(triphenylphosphine) palladium Solvents: Dimethylformamide; 2 h, 110 °C	By: Sandtorv, Alexander H.; et al European Journal of Organic Chemistry (2015), 2015(16), 3506-3512.

Scheme 6 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-1314477 Steps: 1 Yield: 100% 1.1 Reagents: Potassium hydroxide Catalysts: Palladium Solvents: Water; rt 1.2 Reagents: Hydrochloric acid Solvents: Water; rt Experimental Protocols	Dendrimer-Encapsulated Pd Nanoparticles as Aqueous, Room-Temperature Catalysts for the Stille Reaction By: Garcia-Martinez, Joaquin C.; et al Journal of the American Chemical Society (2005), 127(14), 5097-5103.
---	---

Scheme 7 (1 Reaction)


Steps: 1 Yield: 100%



31-176-CAS-6679654 Steps: 1 Yield: 100% 1.1 Catalysts: Tetrakis(triphenylphosphine)palladium Solvents: Dimethylformamide; 12 h, rt	The stereoselective synthesis of (Z)-HfC=CFZnI and stereospecific preparation of (E)-1,2-difluorostyrenes from (Z)-HfC=CFZnI via an unusual Pd(PPh₃)₄-Cu(I)Br co-catalysis approach or (Z)-HfC=CFSnBu₃ By: Liu, Qibo; et al Journal of Fluorine Chemistry (2011), 132(2), 78-87.
---	--

Scheme 8 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-256153 Steps: 1 Yield: 100% 1.1 Reagents: Lithium chloride Catalysts: Dichlorobis(triphenylphosphine)palladium Solvents: Tetrahydrofuran; 480 min, 85 °C Experimental Protocols	Dual Effect of Halides in the Stille Reaction: In Situ Halide Metathesis and Catalyst Stabilization By: Verbeeck, Stefan; et al Chemistry - A European Journal (2010), 16(43), 12831-12837, S12831/1-S12831/17.
--	--

Scheme 9 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-13310233 Steps: 1 Yield: 100% 1.1 Reagents: Potassium hydroxide Catalysts: Palladium, L-Threonyl-L-seryl-L-asparaginy-L-alanyl-L-valyl-L-histidyl-L-prolyl-L-threonyl-L-leucyl-L-arginyl-L-histidyl-L-leucine Solvents: Water; 24 h, 40 °C 1.2 Reagents: Hydrochloric acid Solvents: Water Experimental Protocols	Exploring the mechanism of Stille C-C coupling via peptide-capped Pd nanoparticles results in low temperature reagent selectivity By: Pacardo, Dennis B.; et al Catalysis Science & Technology (2013), 3(3), 745-753.
--	--

Scheme 10 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-392546 Steps: 1 Yield: 100% 1.1 Catalysts: Palladium diacetate Solvents: Dimethylformamide; rt → 100 °C; 25 min, 100 °C; 100 °C → 20 °C	3,4,5-Triarylthiazoles via C-C coupling chemistry By: Christoforou, Irene C.; et al Organic & Biomolecular Chemistry (2007), 5(9), 1381-1390.
--	--

Scheme 11 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-9785363 Steps: 1 Yield: 100% 1.1 Catalysts: Tetrakis(triphenylphosphine)palladium Solvents: Dimethylformamide; 4 h, 100 °C 1.2 Reagents: Potassium fluoride Solvents: Toluene, Water; overnight, rt Experimental Protocols	Palladium-catalyzed reactions of 6-halogeno-1,1'-binaphthyl derivatives. A detailed investigation of structure/reactivity and structure/selectivity relationships By: Feher, Csaba; et al Tetrahedron (2011), 67(34), 6327-6333.
---	---

Scheme 12 (1 Reaction)


Steps: 1 Yield: 100%



31-176-CAS-3779137	Steps: 1 Yield: 100%	Incorporation of phosphole moieties into the side chain of tyrosine and phenylalanine
1.1 Catalysts: Bis(dibenzylideneacetone)palladium Solvents: Tetrahydrofuran; 12 h, 60 °C		By: Bisaro, Fabrice; et al
Experimental Protocols		Synlett (2010), (20), 3081-3085.

Scheme 13 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (59)

Suppliers (38)

31-130-CAS-6867403	Steps: 1 Yield: 100%	Stille Cross-Coupling for the Functionalization of the Imidazole Backbone: Revisit, Improvement, and Applications of the Method
1.1 Catalysts: Cuprous iodide, Tetrakis(triphenylphosphine)palladium Solvents: Dimethylformamide; 2 h, 110 °C		By: Sandtorv, Alexander H.; et al
		European Journal of Organic Chemistry (2015), 2015(16), 3506-3512.

Scheme 14 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (81)

Suppliers (7)

31-130-CAS-5727874	Steps: 1 Yield: 100%	Synthesis of new fused and substituted benzo and pyrido carbazoles via C-2 (het)arylindoles
1.1 Catalysts: Cuprous iodide Solvents: Dimethylacetamide; 30 min, rt		By: Bourderieux, Aurelie; et al
1.2 Catalysts: Tetrakis(triphenylphosphine)palladium; 12 h, 100 °C		Tetrahedron (2008), 64(49), 11012-11019.

Scheme 15 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (51)


Suppliers (95)

Suppliers (115)

31-130-CAS-12983133 Steps: 1 Yield: 100% 1.1 Reagents: Maltose, Dimethylurea, Ammonium chloride Catalysts: Triphenylarsine, Palladium, tris[μ-[(1,2-η:4,5-η)-(1<i>E</i>, 4<i>E</i>)-1,5-diphenyl-1,4-pentadien-3-one]]di-, compd. with trichloromethane (1:1); 6 h, 90 °C; 90 °C → rt	Stille reactions with tetraalkylstannanes and phenyltrialkylstannanes in low melting sugar-urea-salt mixtures By: Imperato, Giovanni; et al Advanced Synthesis & Catalysis (2006), 348(15), 2243-2247.
--	--

Scheme 16 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (63)


Suppliers (51)

Suppliers (41)

31-130-CAS-6921947 Steps: 1 Yield: 100% 1.1 Reagents: Potassium hydroxide Catalysts: Palladium, L-Threonyl-L-seryl-L-asparaginy-L-alanyl-L-valyl-L-histidyl-L-prolyl-L-threonyl-L-leucyl-L-arginy-L-histidyl-L-leucine Solvents: Water; 24 h, 40 °C 1.2 Reagents: Hydrochloric acid Solvents: Water Experimental Protocols	Exploring the mechanism of Stille C-C coupling via peptide-capped Pd nanoparticles results in low temperature reagent selectivity By: Pacardo, Dennis B.; et al Catalysis Science & Technology (2013), 3(3), 745-753.
--	---

Scheme 17 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (87)

Suppliers (4)

31-130-CAS-14860265 Steps: 1 Yield: 100% 1.1 Catalysts: Cuprous iodide Solvents: Dimethylacetamide; 30 min, rt 1.2 Catalysts: Tetrakis(triphenylphosphine)palladium; 6 h, 100 °C	Synthesis of new fused and substituted benzo and pyrido carbazoles via C-2 (het)arylindoles By: Bourderioux, Aurelie; et al Tetrahedron (2008), 64(49), 11012-11019.
---	--

Scheme 18 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (74)


Suppliers (51)

Suppliers (9)

31-130-CAS-2671941 Steps: 1 Yield: 100% 1.1 Reagents: Potassium hydroxide Catalysts: Palladium, L-Threonyl-L-seryl-L-asparaginyl-L-alanyl-L-valyl-L-histidyl-L-prolyl-L-threonyl-L-leucyl-L-arginyl-L-histidyl-L-leucine Solvents: Water; 24 h, 40 °C 1.2 Reagents: Hydrochloric acid Solvents: Water Experimental Protocols	Exploring the mechanism of Stille C-C coupling via peptide-capped Pd nanoparticles results in low temperature reagent selectivity By: Pacardo, Dennis B.; et al Catalysis Science & Technology (2013), 3(3), 745-753.
--	---

Scheme 19 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-2769621 Steps: 1 Yield: 100% 1.1 Catalysts: Palladium diacetate Solvents: Acetonitrile; 100 °C	New regiospecific isothiazole C-C coupling chemistry By: Christoforou, Irene C.; et al Organic & Biomolecular Chemistry (2006), 4(19), 3681-3693.
--	---

Scheme 20 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-9640348 Steps: 1 Yield: 100% 1.1 Reagents: Tetrabutylammonium fluoride Catalysts: Dabco, Palladium diacetate Solvents: 1,4-Dioxane; 16 h, 100 °C Experimental Protocols	Efficient Stille cross-coupling reaction catalyzed by the Pd(OAc)₂/Dabco catalytic system By: Li, Jin-Heng; et al Journal of Organic Chemistry (2005), 70(7), 2832-2834.
---	---

Scheme 21 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-12090905 Steps: 1 Yield: 100%	An efficient Stille cross-coupling reaction catalyzed by Pd(dba)₂/DAB-Cy catalytic system
1.1 Reagents: Potassium fluoride Catalysts: <i>N,N'</i> -1,2-Ethanedithiolenebis[cyclohexanamine], Bis(dibenzylideneacetone)palladium Solvents: 1,4-Dioxane; 16 h, 100 °C	By: Li, Jin-Heng; et al Tetrahedron (2005), 61(30), 7289-7293.

Scheme 22 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-10679561 Steps: 1 Yield: 100%	Nonstoichiometric Stille Coupling Polycondensation for Synthesizing Naphthalene-Diimide-Based π-Conjugated Polymers
1.1 Catalysts: Tri- <i>o</i> -tolylphosphine, Tris(dibenzylideneacetone) dipalladium Solvents: Toluene; 2 h, 90 °C 1.2 Reagents: Hydrochloric acid Solvents: Water	By: Goto, Eisuke; et al ACS Macro Letters (2015), 4(9), 1004-1007.
Experimental Protocols	

Scheme 23 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-7651920 Steps: 1 Yield: 100%	Palladium-catalyzed reactions of 6-halogeno-1,1'-binaphthyl derivatives. A detailed investigation of structure/reactivity and structure/selectivity relationships
1.1 Catalysts: Tetrakis(triphenylphosphine)palladium Solvents: Dimethylformamide; 2 h, 60 °C 1.2 Reagents: Potassium fluoride Solvents: Toluene, Water; overnight, rt	By: Feher, Csaba; et al Tetrahedron (2011), 67(34), 6327-6333.
Experimental Protocols	

Scheme 24 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-9854602	Steps: 1 Yield: 100%	Dendrimer-Encapsulated Pd Nanoparticles as Aqueous, Room-Temperature Catalysts for the Stille Reaction By: Garcia-Martinez, Joaquin C.; et al Journal of the American Chemical Society (2005), 127(14), 5097-5103.
1.1 Reagents: Potassium hydroxide Catalysts: Palladium Solvents: Water; rt		
1.2 Reagents: Hydrochloric acid Solvents: Water; rt		
Experimental Protocols		


Scheme 25 (1 Reaction)

Steps: 1 Yield: 100%




Suppliers (38)

+



Suppliers (66)

→




Suppliers (66)

31-130-CAS-8091448	Steps: 1 Yield: 100%	Dual Effect of Halides in the Stille Reaction: In Situ Halide Metathesis and Catalyst Stabilization By: Verbeeck, Stefan; et al Chemistry - A European Journal (2010), 16(43), 12831-12837, S12831/1-S12831/17.
1.1 Reagents: Lithium chloride Catalysts: Dichlorobis(triphenylphosphine)palladium Solvents: Tetrahydrofuran; 480 min, 85 °C		
Experimental Protocols		


Scheme 26 (4 Reactions)

Steps: 1 Yield: 100%




Suppliers (51)

+



Suppliers (86)

→



Suppliers (88)


31-130-CAS-10534595	Steps: 1 Yield: 100%	Exploring bio-inspired strategies for the production of noble metal nanocatalysts By: Knecht, Marc R.; et al Polymer Preprints (American Chemical Society, Division of Polymer Chemistry) (2010), 51(1), 12-13.
1.1 Reagents: Potassium hydroxide Catalysts: Palladium Solvents: Water; 24 h, 25 °C		
1.2 Reagents: Hydrochloric acid Solvents: Water		
Experimental Protocols		

31-130-CAS-6995716	Steps: 1 Yield: 100%	Biomimetic Synthesis of Pd Nanocatalysts for the Stille Coupling Reaction By: Pacardo, Dennis B.; et al ACS Nano (2009), 3(5), 1288-1296.
1.1 Reagents: Potassium hydroxide Catalysts: Palladium (peptide surface functionalized nanoparticle), L-Threonyl-L-seryl-L-asparaginyl-L-alanyl-L-valyl-L-histidyl-L-prolyl-L-threonyl-L-leucyl-L-arginyl-L-histidyl-L-leucine (surface functionalized palladium nanoparticle with)		
1.2 Reagents: Hydrochloric acid Solvents: Water		
Experimental Protocols		

31-130-CAS-13962126 Steps: 1 Yield: 100% 1.1 Reagents: Sodium bicarbonate, Sodium hydroxide Catalysts: Palladium(1+), diaqua[2-[[bis(1-methylethyl)phosphino-κP]oxy]-4-methylphenyl-κC]-, (SP-4-2)-, sulfate (2:1) Solvents: Water; 720 min, pH 10.5, 100 °C Experimental Protocols	pH-Dependent C-C Coupling Reactions Catalyzed by Water-Soluble Palladacyclic Aqua Catalysts in Water By: Ogo, Seiji; et al Organometallics (2006), 25(2), 331-338.
31-130-CAS-12576854 Steps: 1 Yield: 100% 1.1 Reagents: Potassium hydroxide Catalysts: Palladium Solvents: Water; rt 1.2 Reagents: Hydrochloric acid Solvents: Water; rt Experimental Protocols	Dendrimer-Encapsulated Pd Nanoparticles as Aqueous, Room-Temperature Catalysts for the Stille Reaction By: Garcia-Martinez, Joaquin C.; et al Journal of the American Chemical Society (2005), 127(14), 5097-5103.

Scheme 27 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-3309519 Steps: 1 Yield: 100% 1.1 Catalysts: Tetrakis(triphenylphosphine)palladium Solvents: Dimethylformamide; 2 h, 110 °C	Stille Cross-Coupling for the Functionalization of the Imidazole Backbone: Revisit, Improvement, and Applications of the Method By: Sandtorv, Alexander H.; et al European Journal of Organic Chemistry (2015), 2015(16), 3506-3512.
---	---

Scheme 28 (1 Reaction)


Steps: 1 Yield: 100%



31-130-CAS-11368371 Steps: 1 Yield: 100% 1.1 Catalysts: Dichlorobis(triphenylphosphine)palladium Solvents: Acetonitrile; rt; 45 min, 82 °C Experimental Protocols	Selective Stille Coupling Reactions of 3-Chloro-5-halo (pseudohalo)-4H-1,2,6-thiadiazin-4-ones By: Ioannidou, Heraklidia A.; et al Organic Letters (2011), 13(21), 5886-5889.
--	--

Scheme 29 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (59)

Suppliers (45)

31-130-CAS-15390210

Steps: 1 Yield: 100%

Stille Cross-Coupling for the Functionalization of the Imidazole Backbone: Revisit, Improvement, and Applications of the Method


By: Sandtorv, Alexander H.; et al

European Journal of Organic Chemistry (2015), 2015(16), 3506-3512.

1.1 **Catalysts:** Cuprous iodide, Tetrakis(triphenylphosphine) palladium
Solvents: Dimethylformamide; 2 h, 110 °C

Scheme 30 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (59)

Suppliers (25)

31-130-CAS-4411119

Steps: 1 Yield: 100%

Stille Cross-Coupling for the Functionalization of the Imidazole Backbone: Revisit, Improvement, and Applications of the Method


By: Sandtorv, Alexander H.; et al

European Journal of Organic Chemistry (2015), 2015(16), 3506-3512.

1.1 **Catalysts:** Cuprous iodide, Tetrakis(triphenylphosphine) palladium
Solvents: Dimethylformamide; 2 h, 110 °C

Scheme 31 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (2)

Suppliers (34)

Supplier (1)

31-130-CAS-9259612

Steps: 1 Yield: 100%

Selective Stille Coupling Reactions of 3-Chloro-5-halo (pseudohalo)-4H-1,2,6-thiadiazin-4-ones

By: Ioannidou, Heraklidia A.; et al


Organic Letters (2011), 13(21), 5886-5889.

1.1 **Catalysts:** Dichlorobis(triphenylphosphine)palladium
Solvents: Acetonitrile; rt; 20 min, 82 °C

Experimental Protocols

Scheme 32 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (45)

Suppliers (66)

Suppliers (32)

31-130-CAS-13696432

Steps: 1 Yield: 100%

Dual Effect of Halides in the Stille Reaction: In Situ Halide Metathesis and Catalyst Stabilization

By: Verbeeck, Stefan; et al


Chemistry - A European Journal (2010), 16(43), 12831-12837, S12831/1-S12831/17.

1.1 **Reagents:** Lithium chloride
Catalysts: Dichlorobis(triphenylphosphine)palladium
Solvents: Tetrahydrofuran; 480 min, 85 °C

Experimental Protocols

Scheme 33 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (67)

Suppliers (45)

31-130-CAS-6588726

Steps: 1 Yield: 100%

Synthesis of oligo(thienylfuran)s with thiophene rings at both ends and their structural, electronic, and field-effect properties


By: Miyata, Yasuo; et al

Chemistry - An Asian Journal (2007), 2(12), 1492-1504.

1.1 **Reagents:** Copper oxide (CuO)
Catalysts: Tetrakis(triphenylphosphine)palladium
Solvents: Dimethylformamide; 5 min, 100 °C; 1 h, 100 °C; 100 °C → rt

Scheme 34 (1 Reaction)

Steps: 1 Yield: 100%



Suppliers (41)

Suppliers (81)

Suppliers (116)

31-130-CAS-5462332

Steps: 1 Yield: 100%

Reusable Copper-Catalyzed Cross-Coupling Reactions of Aryl Halides with Organotin in Inexpensive Ionic Liquids

By: Li, Jin-Heng; et al


Journal of Organic Chemistry (2006), 71(19), 7488-7490.

1.1 **Reagents:** Potassium fluoride
Catalysts: Copper oxide (Cu₂O), Tri-*o*-tolylphosphine
Solvents: Tetrabutylammonium bromide; 28 h, 125 - 130 °C

Experimental Protocols


Scheme 35 (2 Reactions)

Steps: 1 Yield: 100%




Suppliers (41)

+



Suppliers (67)

→



Suppliers (82)

31-130-CAS-1616802

Steps: 1 Yield: 100%

- 1.1 **Reagents:** (+)-Lactose, Dimethylurea, Ammonium chloride
Catalysts: Triphenylarsine, Palladium, tris[μ-[(1,2-η:4,5-η)-(1*E*, 4*E*)-1,5-diphenyl-1,4-pentadien-3-one]]di-, compd. with trichloromethane (1:1); 6 h, 90 °C; 90 °C → rt

Stille reactions with tetraalkylstannanes and phenyltrialkylstannanes in low melting sugar-urea-salt mixtures

By: Imperato, Giovanni; et al

Advanced Synthesis & Catalysis (2006), 348(15), 2243-2247.

31-130-CAS-2565944

Steps: 1 Yield: 100%

- 1.1 **Reagents:** Potassium hydroxide
Catalysts: Dabco, Palladium diacetate, Polyethylene glycol
Solvents: Water; 15 h, 80 °C


PEG-400 promoted Pd(OAc)₂/DABCO-catalyzed cross-coupling reactions in aqueous media

By: Li, Jin-Heng; et al

Tetrahedron (2006), 62(1), 31-38.


Scheme 36 (3 Reactions)

Steps: 1 Yield: 100%




Suppliers (41)

+



Suppliers (71)

→



Suppliers (69)

31-614-CAS-38361614

Steps: 1 Yield: 100%

- 1.1 **Reagents:** Cesium fluoride
Catalysts: Bis[2'-(amino-κM)[1,1'-biphenyl]-2-yl-κC]bis[μ-(methanesulfonato-κO:κO)]dipalladium, 3016418-42-8
Solvents: 2-Methyltetrahydrofuran; 6 h, 80 °C

Bulky, electron-rich, renewable: analogues of Beller's phosphine for cross-couplings

By: van der Westhuizen, Danielle; et al

Catalysis Science & Technology (2023), 13(23), 6733-6742.

Experimental Protocols

31-614-CAS-38361659

Steps: 1 Yield: 100%

- 1.1 **Reagents:** Cesium fluoride
Catalysts: Bis[2'-(amino-κM)[1,1'-biphenyl]-2-yl-κC]bis[μ-(methanesulfonato-κO:κO)]dipalladium, Methanesulfonic acid, 1,1,1-trifluoro-, compd. with bis(1,1,5-trimethylheptyl) phosphine (1:1)
Solvents: 2-Methyltetrahydrofuran; 6 h, 80 °C

Bulky, electron-rich, renewable: analogues of Beller's phosphine for cross-couplings

By: van der Westhuizen, Danielle; et al

Catalysis Science & Technology (2023), 13(23), 6733-6742.

Experimental Protocols

31-130-CAS-7705138

Steps: 1 Yield: 100%

- 1.1 **Reagents:** (+)-Lactose, Dimethylurea, Ammonium chloride
Catalysts: Triphenylarsine, Palladium, tris[μ-[(1,2-η:4,5-η)-(1*E*, 4*E*)-1,5-diphenyl-1,4-pentadien-3-one]]di-, compd. with trichloromethane (1:1); 6 h, 90 °C; 90 °C → rt


Stille reactions with tetraalkylstannanes and phenyltrialkylstannanes in low melting sugar-urea-salt mixtures

By: Imperato, Giovanni; et al

Advanced Synthesis & Catalysis (2006), 348(15), 2243-2247.

Scheme 37 (6 Reactions)

Steps: 1 Yield: 100%



Suppliers (51)

Suppliers (94)

Suppliers (101)

31-130-CAS-4627098

Steps: 1 Yield: 100%

Interrogating the catalytic mechanism of nanoparticle mediated Stille coupling reactions employing bio-inspired Pd nanocatalysts

By: Pacardo, Dennis B.; et al

Nanoscale (2011), 3(5), 2194-2201.

- 1.1 **Reagents:** Potassium hydroxide
Catalysts: Palladium, L-Threonyl-L-seryl-L-asparaginy-L-alanyl-L-valyl-L-histidyl-L-prolyl-L-threonyl-L-leucyl-L-arginyl-L-histidyl-L-leucine (palladium nanoparticle surface functionalized with)
Solvents: Water; 24 h, rt

1.2 rt

31-130-CAS-14392741

Steps: 1 Yield: 100%

Effects of the Material Structure on the Catalytic Activity of Peptide-Templated Pd Nanomaterials

By: Bhandari, Rohit; et al

ACS Catalysis (2011), 1(2), 89-98.

- 1.1 **Reagents:** Potassium hydroxide
Solvents: Water; rt
- 1.2 **Reagents:** Oxygen
Catalysts: Palladium; 24 h, rt
- 1.3 **Reagents:** Hydrochloric acid
Solvents: Water; rt

Experimental Protocols

31-130-CAS-9322109

Steps: 1 Yield: 100%

Exploring bio-inspired strategies for the production of noble metal nanocatalysts

By: Knecht, Marc R.; et al

Polymer Preprints (American Chemical Society, Division of Polymer Chemistry) (2010), 51(1), 12-13.

- 1.1 **Reagents:** Potassium hydroxide
Catalysts: Palladium
Solvents: Water; 24 h, 25 °C
- 1.2 **Reagents:** Hydrochloric acid
Solvents: Water

31-130-CAS-2739940

Steps: 1 Yield: 100%

Biomimetic Synthesis of Pd Nanocatalysts for the Stille Coupling Reaction

By: Pacardo, Dennis B.; et al

ACS Nano (2009), 3(5), 1288-1296.

- 1.1 **Reagents:** Potassium hydroxide
Catalysts: Palladium (peptide surface functionalized nanoparticle), L-Threonyl-L-seryl-L-asparaginy-L-alanyl-L-valyl-L-histidyl-L-prolyl-L-threonyl-L-leucyl-L-arginyl-L-histidyl-L-leucine (surface functionalized palladium nanoparticle with)
Solvents: Water; 24 h, rt
- 1.2 **Reagents:** Hydrochloric acid
Solvents: Water

Experimental Protocols

31-130-CAS-1561713

Steps: 1 Yield: 100%

Dendrimer-Encapsulated Pd Nanoparticles versus Palladium Acetate as Catalytic Precursors in the Stille Reaction in Water

By: Bernechea, Maria; et al

Inorganic Chemistry (2009), 48(10), 4491-4496.


- 1.1 **Reagents:** Potassium hydroxide
Catalysts: Palladium, Ethylenediamine-methyl acrylate copolymer
Solvents: Water; 24 h, 40 °C

Experimental Protocols

31-130-CAS-3170646	Steps: 1 Yield: 100%	Dendrimer-Encapsulated Pd Nanoparticles as Aqueous, Room-Temperature Catalysts for the Stille Reaction By: Garcia-Martinez, Joaquin C.; et al Journal of the American Chemical Society (2005), 127(14), 5097-5103.
1.1 Reagents: Potassium hydroxide Catalysts: Palladium Solvents: Water; 15 h, rt		
1.2 Reagents: Hydrochloric acid Solvents: Water; rt		
Experimental Protocols		

Scheme 38 (2 Reactions)

Steps: 1 Yield: 100%



31-130-CAS-6548266	Steps: 1 Yield: 100%	Exploring bio-inspired strategies for the production of noble metal nanocatalysts By: Knecht, Marc R.; et al Polymer Preprints (American Chemical Society, Division of Polymer Chemistry) (2010), 51(1), 12-13.
1.1 Reagents: Potassium hydroxide Catalysts: Palladium Solvents: Water; 24 h, 25 °C		
1.2 Reagents: Hydrochloric acid Solvents: Water		
31-130-CAS-9137194	Steps: 1 Yield: 100%	Biomimetic Synthesis of Pd Nanocatalysts for the Stille Coupling Reaction By: Pacardo, Dennis B.; et al ACS Nano (2009), 3(5), 1288-1296.
1.1 Reagents: Potassium hydroxide Catalysts: Palladium (peptide surface functionalized nanoparticle), L-Threonyl-L-seryl-L-asparaginy-L-alanyl-L-valyl-L-histidyl-L-prolyl-L-threonyl-L-leucyl-L-arginyl-L-histidyl-L-leucine (surface functionalized palladium nanoparticle with)		
1.2 Reagents: Hydrochloric acid Solvents: Water		
Experimental Protocols		