

Initiating Search

November 10, 2025, 11:14 AM

• Search:

CAS Lexicon Search:

Wolff-Kishner reaction - Preferred Concept

Search Tasks

Task	Result Type	View
Returned Reference Results from CAS Lexicon (394)	 References	View Results
Exported: Retrieved Related Reaction Results + Filters (407)	 Reactions	View Results
Filtered By:		
Document Type:	Journal	

Copyright © 2025 American Chemical Society (ACS). All Rights Reserved.


Internal use only. Redistribution is subject to the terms of your CAS SciFinder License Agreement and CAS information Use Policies.

Reactions (50)

[View in CAS SciFinder](#)

Scheme 1 (1 Reaction)

Steps: 1 Yield: 100%



• Na

 Supplier (1)

 Suppliers (38)

31-614-CAS-33997565

Steps: 1 Yield: 100%

Ruthenium(II)-catalyzed deoxygenation of ketones

1.1 **Reagents:** Potassium *tert*-butoxide, Hydrazine hydrate (1:1)
Catalysts: Dimethyl sulfoxide, (*OC*-6-15)-Dichloro[*N*-[2-(diphenylphosphino- κ^P ethyl]-2-pyridinemethanamine- κ^N], κ^N](triphenylphosphine)ruthenium
Solvents: Water; overnight, 75 °C


By: Gui, Ruohua; et al

Chemical Communications (Cambridge, United Kingdom)
 (2022), 58(75), 10572-10575.

Experimental Protocols

Scheme 2 (1 Reaction)

Steps: 1 Yield: 100%



Absolute stereochemistry shown

Absolute stereochemistry shown

 Suppliers (64)

 Suppliers (64)

31-538-CAS-8909179

Steps: 1 Yield: 100%

Optimization of the Wolff-Kishner reaction for the obtainment of tigogenin from hecogenin

1.1 **Reagents:** Hydrazine hydrate (1:1)**Solvents:** Glycol monoethyl ether


By: Moreno, Mayra Reyes; et al

1.2 **Reagents:** Potassium hydroxide

Revista CENIC, Ciencias Quimicas (2001), 32(1), 51-54.

Scheme 3 (2 Reactions)

Steps: 1 Yield: 92-100%



 Suppliers (69)

 Suppliers (280)

31-538-CAS-5254410 Steps: 1 Yield: 100% 1.1 Reagents: Hydrazine Solvents: Glycerol; 5 min, rt → 90 °C 1.2 Reagents: Potassium hydroxide; 5 min, heated	Glycerol as an alternative green medium for carbonyl compound reductions By: Wolfson, Adi; et al Organic Communications (2009), 2(2), 34-41.
31-538-CAS-15006501 Steps: 1 Yield: 92% 1.1 Reagents: Hydrazine hydrate (1:1) Catalysts: Amberlyst A 26 Solvents: Ethylene glycol; 55 min, reflux	An efficient and environmentally friendly Wolff-Kishner reduction catalyzed by strong base anion-exchange resin By: More, Uttam B.; et al Organic Chemistry: An Indian Journal (2008), 4(12), 537-539.

Scheme 4 (1 Reaction)

Steps: 1 Yield: 99%



Double bond geometry shown

Double bond geometry shown


Suppliers (5)

Suppliers (5)

31-614-CAS-33997553 Steps: 1 Yield: 99% 1.1 Reagents: Potassium <i>tert</i> -butoxide, Hydrazine hydrate (1:1) Catalysts: Dimethyl sulfoxide, (OC-6-15)-Dichloro[<i>N</i> -[2-(diphenylphosphino- κP ethyl]-2-pyridinemethanamine- κN^1 , κN^2](triphenylphosphine)ruthenium Solvents: Methanol; overnight, 75 °C	Ruthenium(II)-catalyzed deoxygenation of ketones By: Gui, Ruohua; et al Chemical Communications (Cambridge, United Kingdom) (2022), 58(75), 10572-10575.
Experimental Protocols	

Scheme 5 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (33)

Suppliers (5)

31-614-CAS-33997552 Steps: 1 Yield: 99% 1.1 Reagents: Potassium <i>tert</i> -butoxide, Hydrazine hydrate (1:1) Catalysts: Dimethyl sulfoxide, (OC-6-15)-Dichloro[<i>N</i> -[2-(diphenylphosphino- κP ethyl]-2-pyridinemethanamine- κN^1 , κN^2](triphenylphosphine)ruthenium Solvents: Methanol; overnight, 75 °C	Ruthenium(II)-catalyzed deoxygenation of ketones By: Gui, Ruohua; et al Chemical Communications (Cambridge, United Kingdom) (2022), 58(75), 10572-10575.
Experimental Protocols	

Scheme 6 (1 Reaction)

Steps: 1 Yield: 99%



31-614-CAS-33997548

Steps: 1 Yield: 99%

Ruthenium(II)-catalyzed deoxygenation of ketones

- 1.1 **Reagents:** Potassium *tert*-butoxide, Hydrazine hydrate (1:1)
Catalysts: Dimethyl sulfoxide, (OC-6-15)-Dichloro[*N*-[2-(diphenylphosphino- κ^P ethyl]-2-pyridinemethanamine- κ^N], κ^N](triphenylphosphine)ruthenium
Solvents: Methanol; overnight, 75 °C


By: Gui, Ruohua; et al

Chemical Communications (Cambridge, United Kingdom) (2022), 58(75), 10572-10575.

Experimental Protocols

Scheme 7 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (81)

Suppliers (61)

31-614-CAS-33997551

Steps: 1 Yield: 99%

Ruthenium(II)-catalyzed deoxygenation of ketones

- 1.1 **Reagents:** Potassium *tert*-butoxide, Hydrazine hydrate (1:1)
Catalysts: Dimethyl sulfoxide, (OC-6-15)-Dichloro[*N*-[2-(diphenylphosphino- κ^P ethyl]-2-pyridinemethanamine- κ^N], κ^N](triphenylphosphine)ruthenium
Solvents: Methanol; overnight, 75 °C


By: Gui, Ruohua; et al

Chemical Communications (Cambridge, United Kingdom) (2022), 58(75), 10572-10575.

Experimental Protocols

Scheme 8 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (71)

Suppliers (48)

31-614-CAS-33997563

Steps: 1 Yield: 99%

Ruthenium(II)-catalyzed deoxygenation of ketones

- 1.1 **Reagents:** Potassium *tert*-butoxide, Hydrazine hydrate (1:1)
Catalysts: Dimethyl sulfoxide, (OC-6-15)-Dichloro[*N*-[2-(diphenylphosphino- κ^P ethyl]-2-pyridinemethanamine- κ^N], κ^N](triphenylphosphine)ruthenium
Solvents: Water; overnight, 75 °C


By: Gui, Ruohua; et al

Chemical Communications (Cambridge, United Kingdom) (2022), 58(75), 10572-10575.

Experimental Protocols

Scheme 9 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (80)

Suppliers (32)

31-614-CAS-33997559

Steps: 1 Yield: 99%

Ruthenium(II)-catalyzed deoxygenation of ketones

- 1.1 **Reagents:** Potassium *tert*-butoxide, Hydrazine hydrate (1:1)
Catalysts: Dimethyl sulfoxide, (OC-6-15)-Dichloro[*N*-[2-(diphenylphosphino- κ^P ethyl]-2-pyridinemethanamine- κ^N], κ^N](triphenylphosphine)ruthenium
Solvents: Water; overnight, 75 °C


By: Gui, Ruohua; et al

Chemical Communications (Cambridge, United Kingdom) (2022), 58(75), 10572-10575.

Experimental Protocols

Scheme 10 (1 Reaction)

Steps: 1 Yield: 99%



Suppliers (136)

Suppliers (86)

31-614-CAS-33997549

Steps: 1 Yield: 99%

Ruthenium(II)-catalyzed deoxygenation of ketones

- 1.1 **Reagents:** Potassium *tert*-butoxide, Hydrazine hydrate (1:1)
Catalysts: Dimethyl sulfoxide, (OC-6-15)-Dichloro[*N*-[2-(diphenylphosphino- κ^P ethyl]-2-pyridinemethanamine- κ^N], κ^N](triphenylphosphine)ruthenium
Solvents: Methanol; overnight, 75 °C


By: Gui, Ruohua; et al

Chemical Communications (Cambridge, United Kingdom) (2022), 58(75), 10572-10575.

Experimental Protocols

Scheme 11 (1 Reaction)

Steps: 1 Yield: 98%



Suppliers (35)

Suppliers (39)

31-538-CAS-17259383

Steps: 1 Yield: 98%

From Bis(imidazol-2-yl)methanes to Asymmetrically Substituted Bis(heterocyclo)methanides in Metal Coordination

- 1.1 **Reagents:** Potassium hydroxide, Hydrazine hydrate (1:1)
Solvents: Water; 3 h, 150 °C



By: Dauer, David-R.; et al

European Journal of Inorganic Chemistry (2017), 2017(13), 1966-1978.

Experimental Protocols

Scheme 12 (3 Reactions)

Steps: 1 Yield: 98%



Suppliers (85)

Suppliers (75)

31-538-CAS-14642882

Steps: 1 Yield: 98%

Alternative synthesis of 5-chloromethyl-2,3-dimethoxy-6-methyl-1,4-benzoquinone: a key intermediate for preparing coenzyme Q analogues

By: Wang, Jin; et al

Journal of Chemical Research (2010), 34(12), 724-725.

1.1 **Reagents:** Potassium hydroxide, Hydrazine hydrate (1:1)
Solvents: 2-Methoxyethanol, Water; 2 h, 70 °C; 1 h, 120 °C; 2 h, 150 °C

31-538-CAS-6630263

Steps: 1 Yield: 98%

A green and efficient synthesis of 1-chloromethyl-2,3,4,5-tetramethoxy-6-methylbenzene

By: Wang, Jin; et al

Journal of Chemical Research (2010), 34(12), 717-718.

1.1 **Reagents:** Potassium hydroxide, Hydrazine hydrate (1:1)
Solvents: Ethylene glycol; 2 h, 70 °C; 1 h, 120 °C; 2 h, 150 °C

31-538-CAS-2568319

Steps: 1 Yield: 98%

Asymmetric synthesis of piperazine derivative



By: Sun, Jiaqiang; et al

Xiandai Huagong (2009), 29(11), 59-60.

1.1 **Reagents:** Hydrazine hydrate (1:1)
Solvents: Ethylene glycol, Water; rt → 70 °C
 1.2 **Reagents:** Potassium hydroxide; 3 h, 70 °C; 70 °C → 125 °C; 1.5 h, 125 °C; 125 °C → 150 °C; 3 h, 150 °C

Scheme 13 (1 Reaction)

Steps: 1 Yield: 97%



Suppliers (4)

Suppliers (3)

31-538-CAS-787016

Steps: 1 Yield: 97%

New red-shifted coelenterazine analogues with an extended electronic conjugation

By: Giuliani, Germano; et al



Tetrahedron Letters (2012), 53(38), 5114-5118.

1.1 **Reagents:** Hydrazine hydrate (1:1)
Solvents: Ethylene glycol; 1 h, 100 °C; 100 °C → rt
 1.2 **Reagents:** Potassium hydroxide; 30 min, 240 °C

Experimental Protocols

Scheme 14 (1 Reaction)

Steps: 1 Yield: 97%




Suppliers (84)

Suppliers (80)

31-522-CAS-7399148	Steps: 1 Yield: 97%	One-step double reduction of aryl nitro and carbonyl groups using hydrazine
1.1 Reagents: Hydrazine, Potassium hydroxide; 24 h, 135 °C		By: Diez-Cecilia, Elena; et al
Experimental Protocols		Tetrahedron Letters (2011), 52(50), 6702-6704.

Scheme 15 (1 Reaction)

Steps: 1 Yield: 97%




Suppliers (98)

Suppliers (70)

31-538-CAS-5275280	Steps: 1 Yield: 97%	Microwave induced synthesis of hydrazones and Wolff-Kishner reduction of carbonyl compounds
1.1 Reagents: Hydrazine		By: Gadhwal, Sunil; et al
1.2 Reagents: Potassium hydroxide		Synlett (1999), (10), 1573-1574.
Experimental Protocols		

Scheme 16 (1 Reaction)

Steps: 1 Yield: 96%




Supplier (1)

Suppliers (37)

31-538-CAS-21361448	Steps: 1 Yield: 96%	Synthesis of n-alkylnaphthalenes via semicarbazones
1.1 Reagents: Potassium hydroxide		By: Zengin, Gulay; et al
Solvents: Diethylene glycol; 4.5 h, reflux		Turkish Journal of Chemistry (2006), 30(2), 139-144.

Scheme 17 (1 Reaction)

Steps: 1 Yield: 96%




Suppliers (68)

Suppliers (73)

31-538-CAS-13218023	Steps: 1 Yield: 96%	Rapid Wolff-Kishner reductions in a silicon carbide microreactor
1.1 Reagents: Potassium hydroxide, Hydrazine hydrate (1:1)		By: Newman, Stephen G.; et al
Solvents: Diethylene glycol monoethyl ether; 10 min, 200 psi, 220 °C		Green Chemistry (2014), 16(1), 176-180.
Experimental Protocols		

Scheme 18 (1 Reaction)

Steps: 1 Yield: 95%



Supplier (1)

Suppliers (5)

31-538-CAS-21433158

Steps: 1 Yield: 95%

No Data Available


Chelates of quinolin-8-ol derivatives. II. Synthesis of long-chain n-alkylquinolin-8-ols

By: Uhlemann, Erhard; et al

Zeitschrift fuer Chemie (1983), 23(9), 334.

Scheme 19 (1 Reaction)

Steps: 1 Yield: 95%



Suppliers (4)

Suppliers (4)

31-538-CAS-21432443

Steps: 1 Yield: 95%

No Data Available


Chelates of quinolin-8-ol derivatives. II. Synthesis of long-chain n-alkylquinolin-8-ols

By: Uhlemann, Erhard; et al

Zeitschrift fuer Chemie (1983), 23(9), 334.

Scheme 20 (1 Reaction)

Steps: 1 Yield: 95%



Suppliers (4)

Suppliers (4)

31-538-CAS-21430534

Steps: 1 Yield: 95%

No Data Available


Chelates of quinolin-8-ol derivatives. II. Synthesis of long-chain n-alkylquinolin-8-ols

By: Uhlemann, Erhard; et al

Zeitschrift fuer Chemie (1983), 23(9), 334.

Scheme 21 (1 Reaction)

Steps: 1 Yield: 95%



Suppliers (29)

Suppliers (4)

31-538-CAS-5623854

Steps: 1 Yield: 95%


Simple modification of the Wolff-Kishner reduction

1.1 **Reagents:** Hydrazine hydrate (1:1)
Solvents: Triethylene glycol

By: Huang, Minlon
 Journal of the American Chemical Society (1946), 68, 2487-8.

Scheme 22 (1 Reaction)

Steps: 1 Yield: 95%



Suppliers (97)

Suppliers (64)

31-614-CAS-33997569

Steps: 1 Yield: 95%

Ruthenium(II)-catalyzed deoxygenation of ketones


1.1 **Reagents:** Potassium *tert*-butoxide, Hydrazine hydrate (1:1)
Catalysts: Dimethyl sulfoxide, (OC-6-15)-Dichloro[*N*-[2-(diphenylphosphino- κ^P ethyl]-2-pyridinemethanamine- κ^N], κ^N](triphenylphosphine)ruthenium
Solvents: Methanol; overnight, 75 °C

By: Gui, Ruohua; et al
 Chemical Communications (Cambridge, United Kingdom) (2022), 58(75), 10572-10575.

Experimental Protocols

Scheme 23 (1 Reaction)

Steps: 1 Yield: 95%



Suppliers (34)

Suppliers (8)

31-538-CAS-11660589

Steps: 1 Yield: 95%


The use of a double Wolff-Kishner reduction in the preparation of aceanthrene and aceanthrylene

1.1 **Reagents:** Hydrazine
Solvents: Diethylene glycol

By: Boerrigter, Jos C. Olde; et al
 Recueil des Travaux Chimiques des Pays-Bas (1989), 108(2), 79-80.

Scheme 24 (1 Reaction)

Steps: 1 Yield: 94%



Suppliers (83)

Suppliers (24)

31-614-CAS-33997542

Steps: 1 Yield: 94%

Ruthenium(II)-catalyzed deoxygenation of ketones

1.1 **Reagents:** Potassium *tert*-butoxide, Hydrazine hydrate (1:1)
Catalysts: Dimethyl sulfoxide, (OC-6-15)-Dichloro[*N*-[2-(diphenylphosphino- κ^P ethyl]-2-pyridinemethanamine- κ^N], κ^N](triphenylphosphine)ruthenium
Solvents: Methanol; overnight, 75 °C


By: Gui, Ruohua; et al

Chemical Communications (Cambridge, United Kingdom) (2022), 58(75), 10572-10575.

Experimental Protocols

Scheme 25 (1 Reaction)

Steps: 1 Yield: 94%



Suppliers (33)

Supplier (1)

31-538-CAS-14159829

Steps: 1 Yield: 94%

Comparative study of the reduction of 2- acyl-10-dialkylaminoalkylphenothiazines


1.1 **Reagents:** Hydrazine
Solvents: Ethanol

By: Schmitt, Joseph; et al

Bulletin de la Societe Chimique de France (1961), 340-9.

Scheme 26 (1 Reaction)

Steps: 1 Yield: 94%



Supplier (1)

Suppliers (49)

31-538-CAS-21367123

Steps: 1 Yield: 94%

Synthesis of n-alkylnaphthalenes via semicarbazones


1.1 **Reagents:** Potassium hydroxide
Solvents: Diethylene glycol; 18 h, reflux

By: Zengin, Gulay; et al

Turkish Journal of Chemistry (2006), 30(2), 139-144.

Scheme 27 (1 Reaction)

Steps: 1 Yield: 94%



31-538-CAS-16716915

Steps: 1 Yield: 94%

Synthesis of novel 3,7-diazabicyclo[3.3.1]nonane derivatives


- 1.1 **Reagents:** Hydrazine hydrate (1:1)
Solvents: Triethylene glycol
- 1.2 **Reagents:** Potassium hydroxide; 60 °C; 5 h, 160 - 170 °C; 190 - 200 °C; 1 h, 200 °C → rt
- 1.3 **Reagents:** Water; cooled

By: Malmakova, A. Ye.; et al

Eurasian Chemico-Technological Journal (2014), 16(1), 85-89.

Scheme 28 (1 Reaction)

Steps: 1 Yield: 94%



Suppliers (68)

Suppliers (59)

31-538-CAS-2630198

Steps: 1 Yield: 94%

Rapid Wolff-Kishner reductions in a silicon carbide micror eactor

- 1.1 **Reagents:** Potassium hydroxide, Hydrazine hydrate (1:1)
Solvents: Diethylene glycol monoethyl ether; 3 min, 200 psi, 200 °C


By: Newman, Stephen G.; et al

Green Chemistry (2014), 16(1), 176-180.

Experimental Protocols

Scheme 29 (1 Reaction)

Steps: 1 Yield: 94%



Suppliers (105)

Suppliers (133)

31-538-CAS-7141870

Steps: 1 Yield: 94%

Microwave induced synthesis of hydrazones and Wolff-Kishner reduction of carbonyl compounds

- 1.1 **Reagents:** Hydrazine
- 1.2 **Reagents:** Potassium hydroxide


By: Gadhwal, Sunil; et al

Synlett (1999), (10), 1573-1574.

Experimental Protocols

Scheme 30 (1 Reaction)

Steps: 1 Yield: 93%



Suppliers (4)

Suppliers (79)

31-538-CAS-21226089

Steps: 1 Yield: 93%

Preparation of 2-methyl-1,3-cyclopentanedione


1.1 Solvents: Ethylene glycol

By: Panouse, Jacques J.; et al

Bulletin de la Societe Chimique de France (1955), 1036-9.

Scheme 31 (1 Reaction)

Steps: 1 Yield: 92%



Suppliers (4)

Suppliers (62)

31-538-CAS-9139228

Steps: 1 Yield: 92%


Low-pressure method for Wolff-Kishner reduction
 1.1 Reagents: Sodium, Hydrazine hydrate (1:1)
 Solvents: Diethylene glycol

By: Soffer, Milton D.; et al

Journal of the American Chemical Society (1945), 67, 1435-6.

Scheme 32 (1 Reaction)

Steps: 1 Yield: 92%



31-538-CAS-19883682

Steps: 1 Yield: 92%

Synthetic Studies on Enantioselective Total Synthesis of Cyathane Diterpenoids: Cyrneines A and B, Glaucopine C, and (+)-Allocyathin B₂
 1.1 Reagents: Trifluoroacetic acid, Triethylsilane
 Solvents: Dichloromethane; 0 °C; 2 h, rt

By: Wu, Guo-jie; et al


 1.2 Reagents: Sodium bicarbonate
 Solvents: Water

Journal of Organic Chemistry (2019), 84(6), 3223-3238.


Experimental Protocols

Scheme 33 (3 Reactions)

Steps: 1 Yield: 91-92%



Suppliers (77)



Suppliers (41)

31-538-CAS-759276

Steps: 1 Yield: 92%

**Reduction of steroid ketones and other carbonyl compounds
by modified Wolff-Kishner method**

1.1 Reagents: Hydrazine
Solvents: Diethylene glycol

By: Huang-Minlon

Journal of the American Chemical Society (1949), 71, 3301-3.

31-538-CAS-2875361

Steps: 1 Yield: 91%

**Reduction of steroid ketones and other carbonyl compounds
by modified Wolff-Kishner method**

1.1 Reagents: Hydrazine hydrate (1:1)
Solvents: Diethylene glycol

By: Huang-Minlon

Journal of the American Chemical Society (1949), 71, 3301-3.

31-538-CAS-14157605

Steps: 1 Yield: 91%

**Reduction of steroid ketones and other carbonyl compounds
by modified Wolff-Kishner method**


1.1 Reagents: Hydrazine hydrate (1:1)
Solvents: Diethylene glycol

By: Huang-Minlon


Journal of the American Chemical Society (1949), 71, 3301-3.

Scheme 34 (1 Reaction)

Steps: 1 Yield: 92%



Suppliers (2)



Suppliers (2)

31-538-CAS-9564725

Steps: 1 Yield: 92%

**Efficient Method for the Synthesis of Hetarenoindanones
Based on 3-Arylhetarenes and Their Conversion into
Hetarenoindenes**

1.1 Reagents: Hydrazine
Solvents: Diethylene glycol; 1 h, 80 °C; 1 h, reflux; reflux → rt

By: Kashulin, Igor A.; et al

Journal of Organic Chemistry (2004), 69(16), 5476-5479.

1.2 Reagents: Potassium hydroxide
Solvents: Water; rt; 2 h, reflux

Experimental Protocols

Scheme 35 (1 Reaction)

Steps: 1 Yield: 92%



Suppliers (95)

Suppliers (81)

31-538-CAS-4148621

Steps: 1 Yield: 92%

The simplified Wolff-Kishner reaction with aromatic hydroxy carbonyl compounds


1.1 Reagents: Hydrazine

By: Lock, G.

Monatshefte fuer Chemie (1954), 85, 802-6.

Scheme 36 (1 Reaction)

Steps: 1 Yield: 92%



Suppliers (67)

Suppliers (61)

31-538-CAS-10451801

Steps: 1 Yield: 92%

A new efficient and mild synthesis of 2-oxindoles by one-pot Wolff-Kishner-like reduction of isatin derivatives


1.1 Reagents: Hydrazine hydrate (1:1)

By: Crestini, C.; et al

Synthetic Communications (1994), 24(20), 2835-41.

Scheme 37 (1 Reaction)

Steps: 1 Yield: 92%



Suppliers (88)

Suppliers (46)

31-538-CAS-9547096

Steps: 1 Yield: 92%

Microwave induced synthesis of hydrazones and Wolff-Kishner reduction of carbonyl compounds

1.1 Reagents: Hydrazine

By: Gadhwal, Sunil; et al


1.2 Reagents: Potassium hydroxide

Experimental Protocols

Synlett (1999), (10), 1573-1574.

Scheme 38 (1 Reaction)

Steps: 1 Yield: 92%



Suppliers (101)

Suppliers (86)

31-538-CAS-31834

Steps: 1 Yield: 92%


Modified Wolff-Kishner method

1.1 **Reagents:** Hydrazine hydrate (1:1)
Solvents: Ethylene glycol

By: Huang, Min-Lon; et al
 Hua Hseuh Hseuh Pao (1961), 27(1), 1-9.

Scheme 39 (1 Reaction)

Steps: 1 Yield: 92%



Suppliers (79)

Suppliers (39)

31-538-CAS-7152085

Steps: 1 Yield: 92%


An efficient and environmentally friendly Wolff-Kishner reduction catalyzed by strong base anion-exchange resin

1.1 **Reagents:** Hydrazine hydrate (1:1)
Catalysts: Amberlyst A 26
Solvents: Ethylene glycol; 55 min, reflux

By: More, Uttam B.; et al
 Organic Chemistry: An Indian Journal (2008), 4(12), 537-539.

Scheme 40 (1 Reaction)

Steps: 1 Yield: 91%



Suppliers (4)

Supplier (1)

31-538-CAS-7968309

Steps: 1 Yield: 91%


A New Synthetic Route of 2-Aroyl- and 2-Benzyl-Benzofurans and their Application in the Total Synthesis of a Metabolite Isolated from Dorstenia gigas

1.1 **Reagents:** Hydrazine
Solvents: Ethylene glycol; 15 min, 150 °C
 1.2 **Reagents:** Potassium hydroxide; rt; 4 h, rt → 150 °C
 1.3 **Reagents:** Hydrochloric acid
Solvents: Water; neutralized

By: Correa, Christian; et al
 Australian Journal of Chemistry (2008), 61(12), 991-999.

Scheme 41 (1 Reaction)

Steps: 1 Yield: 91%



Supplier (1)

31-538-CAS-4671377

Steps: 1 Yield: 91%

Synthesis and structure of 7-alkoxyalkyl-3-thia-7-azabicyclo [3.3.1]nonan-9-ones and several of their derivatives

- 1.1 **Reagents:** Potassium hydroxide, Hydrazine hydrate (1:1)
Solvents: Triethylene glycol, Water; 60 °C; 60 °C → 170 °C; 4 h, reflux; reflux → 200 °C; 200 °C → rt


By: Yu, V. K.; et al

Chemistry of Heterocyclic Compounds (New York, NY, United States) (2006), 42(4), 512-519.

Experimental Protocols

Scheme 42 (1 Reaction)

Steps: 1 Yield: 91%



Suppliers (49)

Suppliers (10)

31-538-CAS-2550609

Steps: 1 Yield: 91%

Microwave induced synthesis of hydrazones and Wolff-Kishner reduction of carbonyl compounds

- 1.1 **Reagents:** Hydrazine
 1.2 **Reagents:** Potassium hydroxide


By: Gadhwal, Sunil; et al

Synlett (1999), (10), 1573-1574.

Experimental Protocols

Scheme 43 (1 Reaction)

Steps: 1 Yield: 91%



Suppliers (2)

Suppliers (3)

31-538-CAS-7571102

Steps: 1 Yield: 91%

Facile Multistep Synthesis of Isotru xene and Isotru xenone

- 1.1 **Reagents:** Potassium hydroxide, Hydrazine hydrate (1:1)
Solvents: Ethylene glycol diethyl ether; 24 h, 180 °C; 180 °C → 60 °C

By: Yang, Jye-Shane; et al


Journal of Organic Chemistry (2009), 74(10), 3974-3977.

- 1.2 **Reagents:** Hydrochloric acid
Solvents: Water; 0 °C

Experimental Protocols

Scheme 44 (1 Reaction)

Steps: 1 Yield: 91%



Suppliers (24)

Suppliers (46)

31-399-CAS-4599064

Steps: 1 Yield: 91%

New process for the preparation of 3,5-dihydroxy-1-pentylbenzene

By: Zhang, Long; et al


Research on Chemical Intermediates (2007), 33(6), 535-540.

1.1 **Reagents:** Hydrazine hydrate (1:1); 2 h, reflux; 3 h, reflux;
reflux → rt

Experimental Protocols

Scheme 45 (1 Reaction)

Steps: 1 Yield: 91%



Suppliers (83)

Suppliers (75)

31-538-CAS-12900822

Steps: 1 Yield: 91%

Rapid Wolff-Kishner reductions in a silicon carbide microreactor

By: Newman, Stephen G.; et al

Green Chemistry (2014), 16(1), 176-180.

1.1 **Reagents:** Potassium hydroxide, Hydrazine hydrate (1:1)
Solvents: Diethylene glycol monoethyl ether; 3 min, 200 psi,
200 °C

Experimental Protocols